

LASS2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01385-1

Capture Antibody Information

Catalog Number:
84514-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241898A11
Reactivity:
human
GenBank:
BC010032
Immunogen Catalog Number:
Ag14151

Conjugate:
Unconjugated
Full name:
LAG1 homolog, ceramide synthase 2
Gene ID:
29956

Detection Antibody Information

Catalog Number:
84514-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241898C5
Reactivity:
human
GenBank:
BC010032
Immunogen Catalog Number:
Ag14151

Conjugate:
Unconjugated
Full name:
LAG1 homolog, ceramide synthase 2
Gene ID:
29956

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01385-1 targets LASS2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LASS2 Recombinant antibody, PBS Only (Capture) 84514-2-PBS (241898A11). 100 µg. Concentration 1 mg/ml.

Detection antibody: LASS2 Recombinant antibody, PBS Only (Capture/Detector) 84514-1-PBS (241898C5). 100 µg. Concentration 1 mg/ml.

Alternative LASS2 matched antibody pairs: MP01385-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01385-1, LASS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84514-2-PBS. Detection antibody: 84514-1-PBS. Standard: Ag14151. Range: 0.156-20 ng/mL