

For Research Use Only

MDM4 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01394-1

Capture Antibody Information

Catalog Number:
84534-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241936E4
Reactivity:
human
GenBank:
BC067299
Immunogen Catalog Number:
Ag26678

Conjugate:
Unconjugated
Full name:
Mdm4 p53 binding protein homolog (mouse)
Gene ID:
4194

Detection Antibody Information

Catalog Number:
84534-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241936H5
Reactivity:
human
GenBank:
BC067299
Immunogen Catalog Number:
Ag26678

Conjugate:
Unconjugated
Full name:
Mdm4 p53 binding protein homolog (mouse)
Gene ID:
4194

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01394-1 targets MDM4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MDM4 Recombinant antibody, PBS Only (Capture) 84534-3-PBS (241936E4). 100 μ g. Concentration 1 mg/ml.

Detection antibody: MDM4 Recombinant antibody, PBS Only (Detector) 84534-4-PBS (241936H5). 100 μ g. Concentration 1 mg/ml.

Alternative MDM4 matched antibody pairs: MP01394-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01394-1, MDMX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84534-3-PBS. Detection antibody: 84534-4-PBS. Standard: Ag26678. Range: 0.156-20 ng/mL