

RAD51C Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01399-2**

Capture Antibody Information

Catalog Number: 84537-1-PBS	Clone ID: 241971D4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: RAD51 homolog C (S. cerevisiae)
Isotype: IgG	GenBank: BC107753	Gene ID: 5889
Purification Method: Protein A purification	Immunogen Catalog Number: Ag35070	

Detection Antibody Information

Catalog Number: 84537-3-PBS	Clone ID: 241971F6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: RAD51 homolog C (S. cerevisiae)
Isotype: IgG	GenBank: BC107753	Gene ID: 5889
Purification Method: Protein A purification	Immunogen Catalog Number: Ag35070	

Applications

Tested Applications: Cytometric bead array	Range: 0.469-60 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP01399-2 targets RAD51C in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RAD51C Recombinant antibody, PBS Only (Capture/Detector) 84537-1-PBS (241971D4). 100 µg. Concentration 1 mg/mL.

Detection antibody: RAD51C Recombinant antibody, PBS Only (Detector) 84537-3-PBS (241971F6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01399-2, RAD51C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84537-1-PBS. Detection antibody: 84537-3-PBS. Standard: Ag35070. Range: 0.469-60 ng/mL.