

For Research Use Only

CD99 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01403-1

Capture Antibody Information

Catalog Number:
98174-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241737H1
Reactivity:
human
GenBank:
BC021620

Conjugate:
Unconjugated
Full name:
CD99 molecule
Gene ID:
4267

Detection Antibody Information

Catalog Number:
84580-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241736H8
Reactivity:
human
GenBank:
BC021620

Conjugate:
Unconjugated
Full name:
CD99 molecule
Gene ID:
4267

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01403-1 targets CD99 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Anti-Human CD99 Rabbit Recombinant Antibody, PBS Only 98174-1-PBS (241737H1). 100 µg. Concentration 1 mg/mL.

Detection antibody: CD99 Recombinant antibody, PBS Only (Detector) 84580-1-PBS (241736H8). 100 µg. Concentration 1 mg/mL.

Alternative CD99 matched antibody pairs: MP01403-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01403-1, CD99 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 98174-1-PBS. Detection antibody: 84580-1-PBS. Standard: RP02321. Range: 0.781-100 ng/mL.