

BCAP31 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01418-2

Capture Antibody Information

Catalog Number:
84584-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241978F9
Reactivity:
human
GenBank:
BC014323
Immunogen Catalog Number:
Ag1668

Conjugate:
Unconjugated
Full name:
B-cell receptor-associated protein 31
Gene ID:
10134

Detection Antibody Information

Catalog Number:
84584-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241978H5
Reactivity:
human
GenBank:
BC014323
Immunogen Catalog Number:
Ag1668

Conjugate:
Unconjugated
Full name:
B-cell receptor-associated protein 31
Gene ID:
10134

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01418-2 targets BCAP31 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BCAP31 Recombinant antibody, PBS Only (Capture) 84584-3-PBS (241978F9). 100 µg. Concentration 1 mg/ml.

Detection antibody: BCAP31 Recombinant antibody, PBS Only (Detector) 84584-4-PBS (241978H5). 100 µg. Concentration 1 mg/ml.

Alternative BCAP31 matched antibody pairs: MP01418-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01418-2, BCAP31 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84584-3-PBS. Detection antibody: 84584-4-PBS. Standard: Ag1668. Range: 0.625-80 ng/mL.