

SETD5 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01432-1

Capture Antibody Information

Catalog Number:

84622-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241903B6

Reactivity:

human

GenBank:

BC020956

Immunogen Catalog Number:

Ag28929

Conjugate:

Unconjugated

Full name:

SET domain containing 5

Gene ID:

55209

Detection Antibody Information

Catalog Number:

84622-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241903A3

Reactivity:

human

GenBank:

BC020956

Immunogen Catalog Number:

Ag28929

Conjugate:

Unconjugated

Full name:

SET domain containing 5

Gene ID:

55209

Applications

Tested Applications:

Cytometric bead array

Range:

0.469-60 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01432-1 targets SETD5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SETD5 Recombinant antibody, PBS Only (Capture) 84622-2-PBS (241903B6). 100 µg. Concentration 1 mg/mL.

Detection antibody: SETD5 Recombinant antibody, PBS Only (Detector) 84622-1-PBS (241903A3). 100 µg. Concentration 1 mg/mL.

Alternative SETD5 matched antibody pairs: MP01432-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01432-1, SETD5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84622-2-PBS. Detection antibody: 84622-1-PBS. Standard: Ag28929. Range: 0.469-60 ng/mL