

For Research Use Only

ZNF649 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01434-2

Capture Antibody Information

Catalog Number: 84624-3-PBS	Clone ID: 242036E9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: zinc finger protein 649
Isotype: IgG	GenBank: BC005368	Gene ID: 65251
Purification Method: Protein A purification	Immunogen Catalog Number: Ag22700	

Detection Antibody Information

Catalog Number: 84624-2-PBS	Clone ID: 242036E3	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: zinc finger protein 649
Isotype: IgG	GenBank: BC005368	Gene ID: 65251
Purification Method: Protein A purification	Immunogen Catalog Number: Ag22700	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP01434-2 targets ZNF649 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ZNF649 Recombinant antibody, PBS Only (Capture) 84624-3-PBS (242036E9). 100 µg. Concentration 1 mg/ml.

Detection antibody: ZNF649 Recombinant antibody, PBS Only (Detector) 84624-2-PBS (242036E3). 100 µg. Concentration 1 mg/ml.

Alternative ZNF649 matched antibody pairs: MP01434-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01434-2, ZNF649 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84624-3-PBS. Detection antibody: 84624-2-PBS. Standard: Ag22700. Range: 0.781-100 ng/mL.