

For Research Use Only

# Mouse CD31 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01441-2

## Capture Antibody Information

**Catalog Number:**  
84607-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242041D6  
**Reactivity:**  
mouse  
**GenBank:**  
NM\_001305157.1

**Conjugate:**  
Unconjugated  
**Full name:**  
platelet/endothelial cell adhesion molecule 1  
**Gene ID:**  
18613

## Detection Antibody Information

**Catalog Number:**  
84607-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242041H3  
**Reactivity:**  
mouse  
**GenBank:**  
NM\_001305157.1

**Conjugate:**  
Unconjugated  
**Full name:**  
platelet/endothelial cell adhesion molecule 1  
**Gene ID:**  
18613

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.469-30 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01441-2 targets CD31 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CD31/PECAM-1 Recombinant antibody, PBS Only (Capture) 84607-1-PBS (242041D6). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: Mouse CD31/PECAM-1 Recombinant antibody, PBS Only (Detector) 84607-3-PBS (242041H3). 100  $\mu$ g. Concentration 1 mg/ml.

Alternative CD31 matched antibody pairs: MP01441-1, MP50319-1, MP50319-2, MP50319-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01441-2, MOUSE CD31/PECAM-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84607-1-PBS. Detection antibody: 84607-3-PBS. Standard: Eg0846. Range: 0.469-30 ng/mL