

For Research Use Only

# CD164 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01444-2

## Capture Antibody Information

Catalog Number:  
84598-1-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241950D3

Reactivity:  
human

GenBank:  
NM\_006016.6

Conjugate:  
Unconjugated

Full name:  
CD164 molecule, sialomucin

Gene ID:  
8763

## Detection Antibody Information

Catalog Number:  
84598-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241950G2

Reactivity:  
human

GenBank:  
NM\_006016.6

Conjugate:  
Unconjugated

Full name:  
CD164 molecule, sialomucin

Gene ID:  
8763

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01444-2 targets CD164 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD164 Recombinant antibody, PBS Only (Capture) 84598-1-PBS (241950D3). 100 µg. Concentration 1 mg/ml.

Detection antibody: CD164 Recombinant antibody, PBS Only (Detector) 84598-2-PBS (241950G2). 100 µg. Concentration 1 mg/ml.

Alternative CD164 matched antibody pairs: MP01444-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01444-2, CD164 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84598-1-PBS. Detection antibody: 84598-2-PBS. Standard: Eg1690. Range: 0.781-100 ng/mL.