

INTS6 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01446-2

Capture Antibody Information

Catalog Number:
84610-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242099E6
Reactivity:
human
GenBank:
BC039829
Immunogen Catalog Number:
Ag36146

Conjugate:
Unconjugated
Full name:
integrator complex subunit 6
Gene ID:
26512

Detection Antibody Information

Catalog Number:
84610-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242099D12
Reactivity:
human
GenBank:
BC039829
Immunogen Catalog Number:
Ag36146

Conjugate:
Unconjugated
Full name:
integrator complex subunit 6
Gene ID:
26512

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01446-2 targets INTS6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: INTS6 Recombinant antibody, PBS Only (Capture) 84610-3-PBS (242099E6). 100 µg. Concentration 1 mg/ml.

Detection antibody: INTS6 Recombinant antibody, PBS Only (Detector) 84610-2-PBS (242099D12). 100 µg. Concentration 1 mg/ml.

Alternative INTS6 matched antibody pairs: MP01446-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01446-2, INTS6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84610-3-PBS. Detection antibody: 84610-2-PBS. Standard: Ag36146. Range: 0.313-40 ng/mL