

For Research Use Only

UTX Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01453-2

Capture Antibody Information

Catalog Number: 84614-1-PBS	Clone ID: 242033A9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: ubiquitously transcribed tetrapeptide repeat, X chromosome
Isotype: IgG	GenBank: BC093868	Gene ID: 7403
Purification Method: Protein A purification	Immunogen Catalog Number: Ag36754	

Detection Antibody Information

Catalog Number: 84614-2-PBS	Clone ID: 242033B10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: ubiquitously transcribed tetrapeptide repeat, X chromosome
Isotype: IgG	GenBank: BC093868	Gene ID: 7403
Purification Method: Protein A purification	Immunogen Catalog Number: Ag36754	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP01453-2 targets UTX in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: UTX Recombinant antibody, PBS Only (Capture) 84614-1-PBS (242033A9). 100 µg. Concentration 1 mg/ml.

Detection antibody: UTX Recombinant antibody, PBS Only (Capture/Detector) 84614-2-PBS (242033B10). 100 µg. Concentration 1 mg/ml.

Alternative UTX matched antibody pairs: MP01453-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01453-2, UTX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84614-1-PBS. Detection antibody: 84614-2-PBS. Standard: Ag36754. Range: 0.781-100 ng/mL.