

For Research Use Only

GANC Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01464-1

Capture Antibody Information

Catalog Number:
84676-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242060G6
Reactivity:
human
GenBank:
BC059406
Immunogen Catalog Number:
Ag35067

Conjugate:
Unconjugated
Full name:
glucosidase, alpha; neutral C
Gene ID:
2595

Detection Antibody Information

Catalog Number:
84676-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242060D2
Reactivity:
human
GenBank:
BC059406
Immunogen Catalog Number:
Ag35067

Conjugate:
Unconjugated
Full name:
glucosidase, alpha; neutral C
Gene ID:
2595

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01464-1 targets GANC in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GANC Recombinant antibody, PBS Only (Capture) 84676-2-PBS (242060G6). 100 μ g. Concentration 1 mg/ml.

Detection antibody: GANC Recombinant antibody, PBS Only (Detector) 84676-1-PBS (242060D2). 100 μ g. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01464-1, GANC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84676-2-PBS. Detection antibody: 84676-1-PBS. Standard:null. Range: null.