

# TRIM26 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01471-2

## Capture Antibody Information

**Catalog Number:**  
84661-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242154A8  
**Reactivity:**  
human  
**GenBank:**  
BC024039  
**Immunogen Catalog Number:**  
Ag25546

**Conjugate:**  
Unconjugated  
**Full name:**  
tripartite motif-containing 26  
**Gene ID:**  
7726

## Detection Antibody Information

**Catalog Number:**  
84661-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242154C11  
**Reactivity:**  
human  
**GenBank:**  
BC024039  
**Immunogen Catalog Number:**  
Ag25546

**Conjugate:**  
Unconjugated  
**Full name:**  
tripartite motif-containing 26  
**Gene ID:**  
7726

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01471-2 targets TRIM26 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRIM26 Recombinant antibody, PBS Only (Capture) 84661-1-PBS (242154A8). 100 µg. Concentration 1 mg/ml.

Detection antibody: TRIM26 Recombinant antibody, PBS Only (Detector) 84661-3-PBS (242154C11). 100 µg. Concentration 1 mg/ml.

Alternative TRIM26 matched antibody pairs: MP01471-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP01471-2, TRIM26 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84661-1-PBS. Detection antibody: 84661-3-PBS. Standard: Ag25546. Range: 0.313-40 ng/mL