

RXRB Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01484-1

Capture Antibody Information

Catalog Number:
84656-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242178E1

Reactivity:
human

GenBank:
BC001167

Immunogen Catalog Number:
Ag0263

Conjugate:
Unconjugated

Full name:
retinoid X receptor, beta

Gene ID:
6257

Detection Antibody Information

Catalog Number:
84656-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242178D6

Reactivity:
human

GenBank:
BC001167

Immunogen Catalog Number:
Ag0263

Conjugate:
Unconjugated

Full name:
retinoid X receptor, beta

Gene ID:
6257

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01484-1 targets RXRB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RXRB Recombinant antibody, PBS Only (Capture) 84656-3-PBS (242178E1). 100 µg. Concentration 1 mg/mL.

Detection antibody: RXRB Recombinant antibody, PBS Only (Capture/Detector) 84656-2-PBS (242178D6). 100 µg. Concentration 1 mg/mL.

Alternative RXRB matched antibody pairs: MP01484-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01484-1, RXRB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84656-3-PBS. Detection antibody: 84656-2-PBS. Standard: Ag0263. Range: 0.313-40 ng/mL