

For Research Use Only

RPL11 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01493-1

Capture Antibody Information

Catalog Number:
84670-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242147F7

Reactivity:
human

GenBank:
BC018970

Immunogen Catalog Number:
Ag9356

Conjugate:
Unconjugated

Full name:
ribosomal protein L11

Gene ID:
6135

Detection Antibody Information

Catalog Number:
84670-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242147H5

Reactivity:
human

GenBank:
BC018970

Immunogen Catalog Number:
Ag9356

Conjugate:
Unconjugated

Full name:
ribosomal protein L11

Gene ID:
6135

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01493-1 targets RPL11 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RPL11 Recombinant antibody, PBS Only (Capture) 84670-1-PBS (242147F7). 100 μ g. Concentration 1 mg/ml.

Detection antibody: RPL11 Recombinant antibody, PBS Only (Detector) 84670-2-PBS (242147H5). 100 μ g. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01493-1, RPL11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84670-1-PBS. Detection antibody: 84670-2-PBS. Standard: Ag9356. Range: 0.625-80 ng/mL