

HNF1B Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01499-2

Capture Antibody Information

Catalog Number: 84684-1-PBS	Clone ID: 242120D9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: HNF1 homeobox B
Isotype: IgG	GenBank: BC017714	Gene ID: 6928
Purification Method: Protein A purification	Immunogen Catalog Number: Ag3232	

Detection Antibody Information

Catalog Number: 84684-3-PBS	Clone ID: 242120F8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: HNF1 homeobox B
Isotype: IgG	GenBank: BC017714	Gene ID: 6928
Purification Method: Protein A purification	Immunogen Catalog Number: Ag3232	

Applications

Tested Applications: Cytometric bead array	Range: 1.25-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP01499-2 targets HNF1B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HNF1B Recombinant antibody, PBS Only (Capture) 84684-1-PBS (242120D9). 100 µg. Concentration 1 mg/ml.

Detection antibody: HNF1B Recombinant antibody, PBS Only (Detector) 84684-3-PBS (242120F8). 100 µg. Concentration 1 mg/ml.

Alternative HNF1B matched antibody pairs: MP01499-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01499-2, HNF1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84684-1-PBS. Detection antibody: 84684-3-PBS. Standard: Ag3232. Range: 1.25-40 ng/mL