

For Research Use Only

TAF1D Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01514-2

Capture Antibody Information

Catalog Number:
84712-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242087F4

Reactivity:
human

GenBank:
BC001972

Immunogen Catalog Number:
Ag10014

Conjugate:
Unconjugated

Full name:
TATA box binding protein (TBP)-
associated factor, RNA polymerase I,
D, 41kDa

Gene ID:
79101

Detection Antibody Information

Catalog Number:
84712-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242087C9

Reactivity:
human

GenBank:
BC001972

Immunogen Catalog Number:
Ag10014

Conjugate:
Unconjugated

Full name:
TATA box binding protein (TBP)-
associated factor, RNA polymerase I,
D, 41kDa

Gene ID:
79101

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP01514-2 targets TAF1D in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TAF1D Recombinant antibody, PBS Only (Capture) 84712-2-PBS (242087F4). 100 µg. Concentration 1 mg/ml.

Detection antibody: TAF1D Recombinant antibody, PBS Only (Detector) 84712-1-PBS (242087C9). 100 µg. Concentration 1 mg/ml.

Alternative TAF1D matched antibody pairs: MP01514-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01514-2, TAF1D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84712-2-PBS. Detection antibody: 84712-1-PBS. Standard: Ag10014. Range: 0.625-80 ng/mL