

For Research Use Only

MPHOSPH8 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01529-2

Capture Antibody Information

Catalog Number:

84719-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

242205G2

Reactivity:

human

GenBank:

BC003542

Immunogen Catalog Number:

Ag10268

Conjugate:

Unconjugated

Full name:

M-phase phosphoprotein 8

Gene ID:

54737

Detection Antibody Information

Catalog Number:

84719-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

242205G3

Reactivity:

human

GenBank:

BC003542

Immunogen Catalog Number:

Ag10268

Conjugate:

Unconjugated

Full name:

M-phase phosphoprotein 8

Gene ID:

54737

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01529-2 targets MPHOSPH8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MPHOSPH8 Recombinant antibody, PBS Only (Capture) 84719-1-PBS (242205G2). 100 μ g. Concentration 1 mg/mL.

Detection antibody: MPHOSPH8 Recombinant antibody, PBS Only (Capture/Detector) 84719-2-PBS (242205G3). 100 μ g. Concentration 1 mg/mL.

Alternative MPHOSPH8 matched antibody pairs: MP01529-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01529-2, MPP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84719-1-PBS. Detection antibody: 84719-2-PBS. Standard: Ag10268. Range: 0.625-80 ng/mL