

FOSL2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01531-1

Capture Antibody Information

Catalog Number: 84718-4-PBS	Clone ID: 242225F7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: FOS-like antigen 2
Isotype: IgG	GenBank: BC008899	Gene ID: 2355
Purification Method: Protein A purification	Immunogen Catalog Number: Ag8397	

Detection Antibody Information

Catalog Number: 84718-1-PBS	Clone ID: 242225B1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: FOS-like antigen 2
Isotype: IgG	GenBank: BC008899	Gene ID: 2355
Purification Method: Protein A purification	Immunogen Catalog Number: Ag8397	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP01531-1 targets FOSL2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FOSL2 Recombinant antibody, PBS Only (Capture) 84718-4-PBS (242225F7). 100 µg. Concentration 1 mg/ml.

Detection antibody: FOSL2 Recombinant antibody, PBS Only (Detector) 84718-1-PBS (242225B1). 100 µg. Concentration 1 mg/ml.

Alternative FOSL2 matched antibody pairs: MP01531-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01531-1, FOSL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84718-4-PBS. Detection antibody: 84718-1-PBS. Standard: Ag8397. Range: 0.313-40 ng/mL