

For Research Use Only

# PRDM8 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01532-1

## Capture Antibody Information

Catalog Number:  
84716-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242203D8

Reactivity:  
human

GenBank:  
BC071584

Immunogen Catalog Number:  
Ag23373

Conjugate:  
Unconjugated

Full name:  
PR domain containing 8

Gene ID:  
56978

## Detection Antibody Information

Catalog Number:  
84716-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242203F2

Reactivity:  
human

GenBank:  
BC071584

Immunogen Catalog Number:  
Ag23373

Conjugate:  
Unconjugated

Full name:  
PR domain containing 8

Gene ID:  
56978

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01532-1 targets PRDM8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PRDM8 Recombinant antibody, PBS Only (Capture) 84716-2-PBS (242203D8). 100 µg. Concentration 1 mg/ml.

Detection antibody: PRDM8 Recombinant antibody, PBS Only (Detector) 84716-3-PBS (242203F2). 100 µg. Concentration 1 mg/ml.

Alternative PRDM8 matched antibody pairs: MP01532-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01532-1, PRDM8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84716-2-PBS. Detection antibody: 84716-3-PBS. Standard: Ag23373. Range: 0.781-100 ng/mL.