

TXNDC5 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01536-1

Capture Antibody Information

Catalog Number:
84756-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242195D2

Reactivity:
human

GenBank:
BC001199

Immunogen Catalog Number:
Ag13965

Conjugate:
Unconjugated

Full name:
thioredoxin domain containing 5
(endoplasmic reticulum)

Gene ID:
81567

Detection Antibody Information

Catalog Number:
84756-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242195H8

Reactivity:
human

GenBank:
BC001199

Immunogen Catalog Number:
Ag13965

Conjugate:
Unconjugated

Full name:
thioredoxin domain containing 5
(endoplasmic reticulum)

Gene ID:
81567

Applications

Tested Applications:
Cytometric bead array

Range:
1.56-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01536-1 targets TXNDC5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TXNDC5 Recombinant antibody, PBS Only (Capture) 84756-4-PBS (242195D2). 100 µg. Concentration 1 mg/ml.

Detection antibody: TXNDC5 Recombinant antibody, PBS Only (Detector) 84756-1-PBS (242195H8). 100 µg. Concentration 1 mg/ml.

Alternative TXNDC5 matched antibody pairs: MP01536-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01536-1, TXNDC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84756-4-PBS. Detection antibody: 84756-1-PBS. Standard: Ag13965. Range: 1.56-100 ng/mL.