

KLRB1/CD161 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01539-2

Capture Antibody Information

Catalog Number:
84755-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242140F10

Reactivity:
human

GenBank:
NM_002258.3

Conjugate:
Unconjugated

Full name:
killer cell lectin-like receptor
subfamily B, member 1

Gene ID:
3820

Detection Antibody Information

Catalog Number:
84755-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242140E10

Reactivity:
human

GenBank:
NM_002258.3

Conjugate:
Unconjugated

Full name:
killer cell lectin-like receptor
subfamily B, member 1

Gene ID:
3820

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01539-2 targets KLRB1/CD161 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD161/KLRB1 Recombinant antibody, PBS Only (Capture/Detector) 84755-3-PBS (242140F10). 100 µg. Concentration 1 mg/ml.

Detection antibody: CD161/KLRB1 Recombinant antibody, PBS Only (Detector) 84755-2-PBS (242140E10). 100 µg. Concentration 1 mg/ml.

Alternative KLRB1/CD161 matched antibody pairs: MP01539-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01539-2, CD161/KLRB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84755-3-PBS. Detection antibody: 84755-2-PBS. Standard: Eg1788. Range: 0.156-20 ng/mL.