

For Research Use Only

ACSS2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01551-3

Capture Antibody Information

Catalog Number:
84774-6-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242253H2
Reactivity:
human
GenBank:
BC010141
Immunogen Catalog Number:
Ag9121

Conjugate:
Unconjugated
Full name:
acyl-CoA synthetase short-chain family member 2
Gene ID:
55902

Detection Antibody Information

Catalog Number:
84774-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242253G10
Reactivity:
human
GenBank:
BC010141
Immunogen Catalog Number:
Ag9121

Conjugate:
Unconjugated
Full name:
acyl-CoA synthetase short-chain family member 2
Gene ID:
55902

Applications

Tested Applications:
Sandwich ELISA

Range:
0.313-20 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01551-3 targets ACSS2 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: ACSS2 Recombinant antibody, PBS Only (Capture) 84774-6-PBS (242253H2). 100 μ g. Concentration 1 mg/ml.

Detection antibody: ACSS2 Recombinant antibody, PBS Only (Detector) 84774-5-PBS (242253G10). 100 μ g. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01551-3, Human ACSS2 Recombinant Matched Antibody Pair - PBS only. 84774-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag9121. 84774-5-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



The mean ACSS2 concentration was determined to be 18.54 ng/mL in U-87 MG cell extract based on a 1.00 mg/mL extract load and 23.32 ng/mL in HEK-293 cell extract based on a 3.30 mg/mL extract load.