

MDM2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01554-1

Capture Antibody Information

Catalog Number:
84780-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242080C9
Reactivity:
human
GenBank:
NM_001145337.2
Immunogen Catalog Number:
Ag27255

Conjugate:
Unconjugated
Full name:
Mdm2 p53 binding protein homolog (mouse)
Gene ID:
4193

Detection Antibody Information

Catalog Number:
84780-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242080D3
Reactivity:
human
GenBank:
NM_001145337.2
Immunogen Catalog Number:
Ag27255

Conjugate:
Unconjugated
Full name:
Mdm2 p53 binding protein homolog (mouse)
Gene ID:
4193

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01554-1 targets MDM2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MDM2 Recombinant antibody, PBS Only (Capture) 84780-2-PBS (242080C9). 100 µg. Concentration 1 mg/ml.

Detection antibody: MDM2 Recombinant antibody, PBS Only (Detector) 84780-3-PBS (242080D3). 100 µg. Concentration 1 mg/ml.

Alternative MDM2 matched antibody pairs: MP01554-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01554-1, MDM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84780-2-PBS. Detection antibody: 84780-3-PBS. Standard: Ag27255. Range: 0.313-40 ng/mL