

# PSMD14/POH1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01566-1

## Capture Antibody Information

**Catalog Number:**  
84767-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242160E7  
**Reactivity:**  
human  
**GenBank:**  
BC009524  
**Immunogen Catalog Number:**  
Ag2694

**Conjugate:**  
Unconjugated  
**Full name:**  
proteasome (prosome, macropain)  
26S subunit, non-ATPase, 14  
**Gene ID:**  
10213

## Detection Antibody Information

**Catalog Number:**  
84767-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242160E10  
**Reactivity:**  
human  
**GenBank:**  
BC009524  
**Immunogen Catalog Number:**  
Ag2694

**Conjugate:**  
Unconjugated  
**Full name:**  
proteasome (prosome, macropain)  
26S subunit, non-ATPase, 14  
**Gene ID:**  
10213

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-50 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01566-1 targets PSMD14/POH1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PSMD14/POH1 Recombinant antibody, PBS Only (Capture/Detector) 84767-3-PBS (242160E7). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: PSMD14/POH1 Recombinant antibody, PBS Only (Detector) 84767-2-PBS (242160E10). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP01566-1, PSMD14/POH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84767-3-PBS. Detection antibody: 84767-2-PBS. Standard: Ag2694. Range: 0.781-50 ng/mL.