

PSMD14/POH1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01566-2

Capture Antibody Information

Catalog Number:
84767-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242160E7
Reactivity:
human
GenBank:
BC009524
Immunogen Catalog Number:
Ag2694

Conjugate:
Unconjugated
Full name:
proteasome (prosome, macropain)
26S subunit, non-ATPase, 14
Gene ID:
10213

Detection Antibody Information

Catalog Number:
84767-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242160C12
Reactivity:
human
GenBank:
BC009524
Immunogen Catalog Number:
Ag2694

Conjugate:
Unconjugated
Full name:
proteasome (prosome, macropain)
26S subunit, non-ATPase, 14
Gene ID:
10213

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01566-2 targets PSMD14/POH1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PSMD14/POH1 Recombinant antibody, PBS Only (Capture/Detector) 84767-3-PBS (242160E7). 100 μ g. Concentration 1 mg/mL.

Detection antibody: PSMD14/POH1 Recombinant antibody, PBS Only (Capture/Detector) 84767-1-PBS (242160C12). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01566-2, PSMD14/POH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84767-3-PBS. Detection antibody: 84767-1-PBS. Standard: Ag2694. Range: 0.781-50 ng/mL.