

# MCM3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01593-3

## Capture Antibody Information

**Catalog Number:**  
84792-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242103G1  
**Reactivity:**  
human  
**GenBank:**  
BC001626  
**Immunogen Catalog Number:**  
Ag7966

**Conjugate:**  
Unconjugated  
**Full name:**  
minichromosome maintenance complex component 3  
**Gene ID:**  
4172

## Detection Antibody Information

**Catalog Number:**  
84792-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242103A9  
**Reactivity:**  
human  
**GenBank:**  
BC001626  
**Immunogen Catalog Number:**  
Ag7966

**Conjugate:**  
Unconjugated  
**Full name:**  
minichromosome maintenance complex component 3  
**Gene ID:**  
4172

## Applications

**Tested Applications:**  
Sandwich ELISA

**Range:**  
0.391-25 ng/mL (Sandwich ELISA)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01593-3 targets MCM3 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: MCM3 Recombinant antibody, PBS Only (Capture/Detector) 84792-3-PBS (242103G1). 100 µg. Concentration 1 mg/ml.

Detection antibody: MCM3 Recombinant antibody, PBS Only (Detector) 84792-4-PBS (242103A9). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Sandwich ELISA standard curve of MP01593-3, Human MCM3 Recombinant Matched Antibody Pair - PBS only. 84792-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag7966. 84792-4-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL



The mean MCM3 concentration was determined to be 45.31 ng/mL in HepG2 cell extract based on a 1.00 mg/mL extract load and 51.56 ng/mL in L02 cell extract based on a 1.30 mg/mL extract load.