

EMR1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01597-2

Capture Antibody Information

Catalog Number:
84786-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241999B5
Reactivity:
human
GenBank:
NM_001974.5

Conjugate:
Unconjugated
Full name:
egf-like module containing, mucin-like, hormone receptor-like 1
Gene ID:
2015

Detection Antibody Information

Catalog Number:
84786-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241999C1
Reactivity:
human
GenBank:
NM_001974.5

Conjugate:
Unconjugated
Full name:
egf-like module containing, mucin-like, hormone receptor-like 1
Gene ID:
2015

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01597-2 targets EMR1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EMR1/ADGRE1/(F4/80) Recombinant antibody, PBS Only (Capture/Detector) 84786-1-PBS (241999B5). 100 µg. Concentration 1 mg/ml.

Detection antibody: EMR1/ADGRE1/(F4/80) Recombinant antibody, PBS Only (Detector) 84786-2-PBS (241999C1). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01597-2, EMR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84786-1-PBS. Detection antibody: 84786-2-PBS. Standard: Eg1716. Range: 0.78-100 ng/mL