

For Research Use Only

# DNMT3A Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01598-1

## Capture Antibody Information

Catalog Number:  
81474-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242257E3

Reactivity:  
human

GenBank:  
BC043617

Immunogen Catalog Number:  
Ag15111

Conjugate:  
Unconjugated

Full name:  
DNA (cytosine-5-)-methyltransferase  
3 alpha

Gene ID:  
1788

## Detection Antibody Information

Catalog Number:  
81474-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242257D3

Reactivity:  
human

GenBank:  
BC043617

Immunogen Catalog Number:  
Ag15111

Conjugate:  
Unconjugated

Full name:  
DNA (cytosine-5-)-methyltransferase  
3 alpha

Gene ID:  
1788

## Applications

Tested Applications:  
Cytometric bead array

Range:  
2.5-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01598-1 targets DNMT3A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DNMT3A Recombinant antibody, PBS Only (Capture) 81474-3-PBS (242257E3). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: DNMT3A Recombinant antibody, PBS Only (Capture/Detector) 81474-2-PBS (242257D3). 100  $\mu$ g. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01598-1, DNMT3A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81474-3-PBS. Detection antibody: 81474-2-PBS. Standard: Ag15111. Range: 2.5-80 ng/mL.