

For Research Use Only

# AMACR/p504S Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01623-2

## Capture Antibody Information

**Catalog Number:**  
84891-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242456D1  
**Reactivity:**  
human  
**GenBank:**  
BC009471  
**Immunogen Catalog Number:**  
Ag8720

**Conjugate:**  
Unconjugated  
**Full name:**  
alpha-methylacyl-CoA racemase  
**Gene ID:**  
23600

## Detection Antibody Information

**Catalog Number:**  
84891-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242456G6  
**Reactivity:**  
human  
**GenBank:**  
BC009471  
**Immunogen Catalog Number:**  
Ag8720

**Conjugate:**  
Unconjugated  
**Full name:**  
alpha-methylacyl-CoA racemase  
**Gene ID:**  
23600

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.563-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01623-2 targets AMACR/p504S in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AMACR/p504S Recombinant antibody, PBS Only (Capture) 84891-2-PBS (242456D1). 100 µg. Concentration 1 mg/ml.

Detection antibody: AMACR/p504S Recombinant antibody, PBS Only (Detector) 84891-3-PBS (242456G6). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01623-2, AMACR/p504S Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84891-2-PBS. Detection antibody: 84891-3-PBS. Standard: Ag8720. Range: 1.563-100 ng/mL.