

For Research Use Only

Integrin Alpha V Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01650-2

Capture Antibody Information

Catalog Number:
84883-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242376C5
Reactivity:
human
GenBank:
BC136442
Immunogen Catalog Number:
Ag25825

Conjugate:
Unconjugated
Full name:
integrin, alpha V (vitronectin receptor,
alpha polypeptide, antigen CD51)
Gene ID:
3685

Detection Antibody Information

Catalog Number:
84883-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242376D11
Reactivity:
human
GenBank:
BC136442
Immunogen Catalog Number:
Ag25825

Conjugate:
Unconjugated
Full name:
integrin, alpha V (vitronectin receptor,
alpha polypeptide, antigen CD51)
Gene ID:
3685

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP01650-2 targets Integrin Alpha V in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ITGAV Recombinant antibody, PBS Only (Capture) 84883-1-PBS (242376C5). 100 µg. Concentration 1 mg/ml.

Detection antibody: ITGAV Recombinant antibody, PBS Only (Detector) 84883-2-PBS (242376D11). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01650-2, Integrin Alpha V Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84883-1-PBS. Detection antibody: 84883-2-PBS. Standard: Ag25825. Range: 0.781-100 ng/mL.