

For Research Use Only

# NDUFS2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01653-1

## Capture Antibody Information

Catalog Number:  
84880-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242238H6

Reactivity:  
human

GenBank:  
BC000170

Immunogen Catalog Number:  
Ag28008

Conjugate:  
Unconjugated

Full name:  
NADH dehydrogenase (ubiquinone)  
Fe-S protein 2, 49kDa (NADH-  
coenzyme Q reductase)

Gene ID:  
4720

## Detection Antibody Information

Catalog Number:  
84880-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242238H7

Reactivity:  
human

GenBank:  
BC000170

Immunogen Catalog Number:  
Ag28008

Conjugate:  
Unconjugated

Full name:  
NADH dehydrogenase (ubiquinone)  
Fe-S protein 2, 49kDa (NADH-  
coenzyme Q reductase)

Gene ID:  
4720

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.938-60 ng/mL (Cytometric Bead  
Array)

Recommended Dilutions:  
It is recommended that this reagent  
should be titrated in each testing  
system to obtain optimal results.

## Product Information

MP01653-1 targets NDUFS2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NDUFS2 Recombinant antibody, PBS Only (Capture) 84880-2-PBS (242238H6). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: NDUFS2 Recombinant antibody, PBS Only (Capture/Detector) 84880-3-PBS (242238H7). 100  $\mu$ g. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01653-1, NDUFS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84880-2-PBS. Detection antibody: 84880-3-PBS. Standard: Ag28008. Range: 0.938-60 ng/mL.