

For Research Use Only

# BCL3 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01655-2

## Capture Antibody Information

Catalog Number:  
84882-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242233G12

Reactivity:  
human

GenBank:  
BC064993

Immunogen Catalog Number:  
Ag13375

Conjugate:  
Unconjugated

Full name:  
B-cell CLL/lymphoma 3

Gene ID:  
602

## Detection Antibody Information

Catalog Number:  
84882-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242233E1

Reactivity:  
human

GenBank:  
BC064993

Immunogen Catalog Number:  
Ag13375

Conjugate:  
Unconjugated

Full name:  
B-cell CLL/lymphoma 3

Gene ID:  
602

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01655-2 targets BCL3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BCL3 Recombinant antibody, PBS Only (Capture/Detector) 84882-3-PBS (242233G12). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: BCL3 Recombinant antibody, PBS Only (Detector) 84882-2-PBS (242233E1). 100  $\mu$ g. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01655-2, BCL3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84882-3-PBS. Detection antibody: 84882-2-PBS. Standard: Ag13375. Range: 1.563-100 ng/mL.