

ARK5 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01658-1

Capture Antibody Information

Catalog Number:
84893-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242411F7

Reactivity:
human

GenBank:
BC152462

Immunogen Catalog Number:
Ag18832

Conjugate:
Unconjugated

Full name:
NUAK family, SNF1-like kinase, 1

Gene ID:
9891

Detection Antibody Information

Catalog Number:
84893-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242411A2

Reactivity:
human

GenBank:
BC152462

Immunogen Catalog Number:
Ag18832

Conjugate:
Unconjugated

Full name:
NUAK family, SNF1-like kinase, 1

Gene ID:
9891

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA

Range:
0.781-100 ng/mL (Cytometric Bead Array)
0.156-10 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01658-1 targets ARK5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array, Sandwich ELISA.

Capture antibody: NUAK1 Recombinant antibody, PBS Only (Capture) 84893-4-PBS (242411F7). 100 µg. Concentration 1 mg/ml.

Detection antibody: NUAK1 Recombinant antibody, PBS Only (Detector) 84893-1-PBS (242411A2). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.


The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only


Selected Validation Data



Sandwich ELISA standard curve of MP01658-1, Human ARK5 Recombinant Matched Antibody Pair - PBS only. 84893-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag18832. 84893-1-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



The mean ARK5 concentration was determined to be 145.81 ng/mL in PC-3 cell extract based on a 1.50 mg/mL extract load.



Cytometric bead array standard curve of MP01658-1, ARK5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84893-4-PBS. Detection antibody: 84893-1-PBS. Standard: Ag18832. Range: 0.781-100 ng/mL