

MICA Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01661-2

Capture Antibody Information

Catalog Number:

84850-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

242004C9

Reactivity:

human

GenBank:

BC016929

Conjugate:

Unconjugated

Full name:

MHC class I polypeptide-related sequence A

Gene ID:

100507436

Detection Antibody Information

Catalog Number:

84850-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

242004C11

Reactivity:

human

GenBank:

BC016929

Conjugate:

Unconjugated

Full name:

MHC class I polypeptide-related sequence A

Gene ID:

100507436

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01661-2 targets MICA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MICA Recombinant antibody, PBS Only (Capture/Detector) 84850-2-PBS (242004C9). 100 µg. Concentration 1 mg/ml.

Detection antibody: MICA Recombinant antibody, PBS Only (Detector) 84850-1-PBS (242004C11). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01661-2, MICA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84850-2-PBS. Detection antibody: 84850-1-PBS. Standard: Eg0181. Range: 0.625-80 ng/mL