

# Mouse CD25/IL-2RA Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01669-1**

## Capture Antibody Information

<b>Catalog Number:</b> 84925-3-PBS	<b>Clone ID:</b> 242494G8	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> mouse	<b>Full name:</b> interleukin 2 receptor, alpha chain
<b>Isotype:</b> IgG	<b>GenBank:</b> NM_008367.3	<b>Gene ID:</b> 16184
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Eg3087	

## Detection Antibody Information

<b>Catalog Number:</b> 84925-1-PBS	<b>Clone ID:</b> 242494A6	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> mouse	<b>Full name:</b> interleukin 2 receptor, alpha chain
<b>Isotype:</b> IgG	<b>GenBank:</b> NM_008367.3	<b>Gene ID:</b> 16184
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Eg3087	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.625-80 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

## Product Information

MP01669-1 targets CD25/IL-2RA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CD25/IL-2RA Recombinant antibody, PBS Only (Capture) 84925-3-PBS (242494G8). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse CD25/IL-2RA Recombinant antibody, PBS Only (Detector) 84925-1-PBS (242494A6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP01669-1, MOUSE CD25/IL2RA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84925-3-PBS. Detection antibody: 84925-1-PBS. Standard: Eg3087. Range: 0.625-80 ng/mL