

Mouse CD25/IL-2RA Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01669-2**

Capture Antibody Information

Catalog Number: 84925-1-PBS	Clone ID: 242494A6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: mouse	Full name: interleukin 2 receptor, alpha chain
Isotype: IgG	GenBank: NM_008367.3	Gene ID: 16184
Purification Method: Protein A purification	Immunogen Catalog Number: Eg3087	

Detection Antibody Information

Catalog Number: 84925-2-PBS	Clone ID: 242494F8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: mouse	Full name: interleukin 2 receptor, alpha chain
Isotype: IgG	GenBank: NM_008367.3	Gene ID: 16184
Purification Method: Protein A purification	Immunogen Catalog Number: Eg3087	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-80 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP01669-2 targets CD25/IL-2RA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CD25/IL-2RA Recombinant antibody, PBS Only (Detector) 84925-1-PBS (242494A6). 100 μ g. Concentration 1 mg/mL.

Detection antibody: Mouse CD25/IL-2RA Recombinant antibody, PBS Only 84925-2-PBS (242494F8). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01669-2, MOUSE CD25/IL2RA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84925-1-PBS. Detection antibody: 84925-2-PBS. Standard: Eg3087. Range: 0.625-80 ng/mL