

TEAD4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01682-2

Capture Antibody Information

Catalog Number:
84932-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242496B7

Reactivity:
human

GenBank:
BC015497

Immunogen Catalog Number:
Ag3059

Conjugate:
Unconjugated

Full name:
TEA domain family member 4

Gene ID:
7004

Detection Antibody Information

Catalog Number:
84932-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242496E1

Reactivity:
human

GenBank:
BC015497

Immunogen Catalog Number:
Ag3059

Conjugate:
Unconjugated

Full name:
TEA domain family member 4

Gene ID:
7004

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01682-2 targets TEAD4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TEAD4 Recombinant antibody, PBS Only (Capture/Detector) 84932-2-PBS (242496B7). 100 µg. Concentration 1 mg/ml.

Detection antibody: TEAD4 Recombinant antibody, PBS Only (Detector) 84932-3-PBS (242496E1). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01682-2, TEAD4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84932-2-PBS. Detection antibody: 84932-3-PBS. Standard: Ag3059. Range: 0.781-100 ng/mL.