

For Research Use Only

IL-5RA/CD125 Recombinant Matched Antibody Pair, PBS Only



Catalog Number: MP01699-3

Capture Antibody Information

Catalog Number:
84948-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242234B1

Reactivity:
human

GenBank:
NM_175726.4

Conjugate:
Unconjugated

Full name:
interleukin 5 receptor, alpha

Gene ID:
3568

Detection Antibody Information

Catalog Number:
84948-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242234H3

Reactivity:
human

GenBank:
NM_175726.4

Conjugate:
Unconjugated

Full name:
interleukin 5 receptor, alpha

Gene ID:
3568

Applications

Tested Applications:
Sandwich ELISA

Range:
0.391-25 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01699-3 targets IL-5RA/CD125 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: IL-5 R alpha/CD125 Recombinant antibody, PBS Only (Capture) 84948-1-PBS (242234B1). 100 μ g. Concentration 1 mg/ml.

Detection antibody: IL-5 R alpha/CD125 Recombinant antibody, PBS Only (Detector) 84948-4-PBS (242234H3). 100 μ g. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP01699-3, Human IL-5RA/CD125 Recombinant Matched Antibody Pair - PBS only. 84948-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg1665. 84948-4-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL



Serum of six individual healthy human donors was measured. The IL-5RA/CD125 concentration of detected samples was determined to be 1.36 ng/mL with a range of 0.79-1.87 ng/mL.



Plasma of six individual healthy human donors was measured. The IL-5RA/CD125 concentration of detected samples was determined to be 1.86 ng/mL with a range of 0.39-6.09 ng/mL.