

For Research Use Only

HEXIM1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01740-2

Capture Antibody Information

Catalog Number:
85028-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242587F1
Reactivity:
human
GenBank:
BC006460
Immunogen Catalog Number:
Ag8144

Conjugate:
Unconjugated
Full name:
hexamethylene bis-acetamide inducible 1
Gene ID:
10614

Detection Antibody Information

Catalog Number:
85028-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242587A2
Reactivity:
human
GenBank:
BC006460
Immunogen Catalog Number:
Ag8144

Conjugate:
Unconjugated
Full name:
hexamethylene bis-acetamide inducible 1
Gene ID:
10614

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01740-2 targets HEXIM1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HEXIM1 Recombinant antibody, PBS Only (Capture) 85028-3-PBS (242587F1). 100 µg. Concentration 1 mg/mL.

Detection antibody: HEXIM1 Recombinant antibody, PBS Only (Detector) 85028-1-PBS (242587A2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01740-2, HEXIM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85028-3-PBS. Detection antibody: 85028-1-PBS. Standard: Ag8144. Range: 0.156-20 ng/mL