

For Research Use Only

# HEXIM1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01740-3

## Capture Antibody Information

**Catalog Number:**  
85028-5-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242587D9  
**Reactivity:**  
human  
**GenBank:**  
BC006460  
**Immunogen Catalog Number:**  
Ag8144

**Conjugate:**  
Unconjugated  
**Full name:**  
hexamethylene bis-acetamide inducible 1  
**Gene ID:**  
10614

## Detection Antibody Information

**Catalog Number:**  
85028-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242587D5  
**Reactivity:**  
human  
**GenBank:**  
BC006460  
**Immunogen Catalog Number:**  
Ag8144

**Conjugate:**  
Unconjugated  
**Full name:**  
hexamethylene bis-acetamide inducible 1  
**Gene ID:**  
10614

## Applications

**Tested Applications:**  
Sandwich ELISA

**Range:**  
0.391-25 ng/mL (Sandwich ELISA)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01740-3 targets HEXIM1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: HEXIM1 Recombinant antibody, PBS Only (Capture) 85028-5-PBS (242587D9). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: HEXIM1 Recombinant antibody, PBS Only (Detector) 85028-4-PBS (242587D5). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP01740-3, Human HEXIM1 Recombinant Matched Antibody Pair - PBS only. 85028-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag8144. 85028-4-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL



The mean HEXIM1 concentration was determined to be 122.36 ng/mL in HeLa cell extract based on a 1.50 mg/mL extract load, 285.55 ng/mL in MCF-7 cell extract based on a 3.50 mg/mL extract load and 310.04 ng/mL in Jurkat cell extract based on a 3.10 mg/mL extract load.