

For Research Use Only

MLYCD Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01785-1

Capture Antibody Information

Catalog Number:
85080-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242753H2

Reactivity:
human

GenBank:
BC000286

Immunogen Catalog Number:
Ag7334

Conjugate:
Unconjugated

Full name:
malonyl-CoA decarboxylase

Gene ID:
23417

Detection Antibody Information

Catalog Number:
85080-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242753E4

Reactivity:
human

GenBank:
BC000286

Immunogen Catalog Number:
Ag7334

Conjugate:
Unconjugated

Full name:
malonyl-CoA decarboxylase

Gene ID:
23417

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01785-1 targets MLYCD in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MLYCD Recombinant antibody, PBS Only (Capture) 85080-3-PBS (242753H2). 100 μ g. Concentration 1 mg/ml.

Detection antibody: MLYCD Recombinant antibody, PBS Only (Detector) 85080-2-PBS (242753E4). 100 μ g. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01785-1, MLYCD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85080-3-PBS. Detection antibody: 85080-2-PBS. Standard: Ag7334. Range: 0.781-100 ng/mL.