

ASIC1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01789-1

Capture Antibody Information

Catalog Number:
85089-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242749B11
Reactivity:
human
GenBank:
NM_001095
Immunogen Catalog Number:
Ag25941

Conjugate:
Unconjugated
Full name:
amiloride-sensitive cation channel 2, neuronal
Gene ID:
41

Detection Antibody Information

Catalog Number:
85089-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242749H5
Reactivity:
human
GenBank:
NM_001095
Immunogen Catalog Number:
Ag25941

Conjugate:
Unconjugated
Full name:
amiloride-sensitive cation channel 2, neuronal
Gene ID:
41

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01789-1 targets ASIC1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ASIC1 Recombinant antibody, PBS Only (Capture) 85089-1-PBS (242749B11). 100 µg. Concentration 1 mg/ml.

Detection antibody: ASIC1 Recombinant antibody, PBS Only (Detector) 85089-3-PBS (242749H5). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01789-1, ASIC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85089-1-PBS. Detection antibody: 85089-3-PBS. Standard: Ag25941. Range: 0.313-40 ng/mL