

For Research Use Only

DDX46 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01807-2

Capture Antibody Information

Catalog Number:
85065-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242715E10
Reactivity:
human
GenBank:
BC012304
Immunogen Catalog Number:
Ag10342

Conjugate:
Unconjugated
Full name:
DEAD (Asp-Glu-Ala-Asp) box polypeptide 46
Gene ID:
9879

Detection Antibody Information

Catalog Number:
85065-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242715E1
Reactivity:
human
GenBank:
BC012304
Immunogen Catalog Number:
Ag10342

Conjugate:
Unconjugated
Full name:
DEAD (Asp-Glu-Ala-Asp) box polypeptide 46
Gene ID:
9879

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01807-2 targets DDX46 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DDX46 Recombinant antibody, PBS Only (Capture) 85065-3-PBS (242715E10). 100 μ g. Concentration 1 mg/ml.

Detection antibody: DDX46 Recombinant antibody, PBS Only (Detector) 85065-2-PBS (242715E1). 100 μ g. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01807-2, DDX46 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85065-3-PBS. Detection antibody: 85065-2-PBS. Standard: Ag10342. Range: 0.781-100 ng/mL.