

For Research Use Only

CD27 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01811-2

Capture Antibody Information

Catalog Number:
85066-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242731G8

Reactivity:
human

GenBank:
BC012160

Conjugate:
Unconjugated

Full name:
CD27 molecule

Gene ID:
939

Detection Antibody Information

Catalog Number:
85066-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242731B7

Reactivity:
human

GenBank:
BC012160

Conjugate:
Unconjugated

Full name:
CD27 molecule

Gene ID:
939

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01811-2 targets CD27 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD27 Recombinant antibody, PBS Only (Capture) 85066-3-PBS (242731G8). 100 µg. Concentration 1 mg/mL.

Detection antibody: CD27 Recombinant antibody, PBS Only (Detector) 85066-1-PBS (242731B7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01811-2, CD27 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85066-3-PBS. Detection antibody: 85066-1-PBS. Standard: Eg1644. Range: 0.781-100 ng/mL.