

For Research Use Only

# STAT2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01853-1

## Capture Antibody Information

Catalog Number:  
80061-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242685E1

Reactivity:  
human

GenBank:  
BC051284

Immunogen Catalog Number:  
Ag10168

Conjugate:  
Unconjugated

Full name:  
signal transducer and activator of transcription 2, 113kDa

Gene ID:  
6773

## Detection Antibody Information

Catalog Number:  
80061-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
242685E6

Reactivity:  
human

GenBank:  
BC051284

Immunogen Catalog Number:  
Ag10168

Conjugate:  
Unconjugated

Full name:  
signal transducer and activator of transcription 2, 113kDa

Gene ID:  
6773

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01853-1 targets STAT2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: STAT2 Recombinant antibody, PBS Only (Capture) 80061-3-PBS (242685E1). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: STAT2 Recombinant antibody, PBS Only (Detector) 80061-4-PBS (242685E6). 100  $\mu$ g. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01853-1, STAT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80061-3-PBS. Detection antibody: 80061-4-PBS. Standard: Ag10168. Range: 1.563-200 ng/mL.