

For Research Use Only

ETS2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01858-2

Capture Antibody Information

Catalog Number:
85113-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242817H3
Reactivity:
human
GenBank:
BC017040
Immunogen Catalog Number:
Ag2929

Conjugate:
Unconjugated
Full name:
v-ets erythroblastosis virus E26 oncogene homolog 2 (avian)
Gene ID:
2114

Detection Antibody Information

Catalog Number:
85113-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242817A10
Reactivity:
human
GenBank:
BC017040
Immunogen Catalog Number:
Ag2929

Conjugate:
Unconjugated
Full name:
v-ets erythroblastosis virus E26 oncogene homolog 2 (avian)
Gene ID:
2114

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01858-2 targets ETS2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ETS2 Recombinant antibody, PBS Only (Capture) 85113-3-PBS (242817H3). 100 μ g. Concentration 1 mg/mL.

Detection antibody: ETS2 Recombinant antibody, PBS Only (Detector) 85113-1-PBS (242817A10). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01858-2, ETS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85113-3-PBS. Detection antibody: 85113-1-PBS. Standard: Ag2929. Range: 1.563-200 ng/mL.