

For Research Use Only

PCBP3 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01894-1

Capture Antibody Information

Catalog Number:
85180-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242791E7

Reactivity:
human

GenBank:
BC012061

Immunogen Catalog Number:
Ag9202

Conjugate:
Unconjugated

Full name:
poly(rC) binding protein 3

Gene ID:
54039

Detection Antibody Information

Catalog Number:
85180-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
242791F11

Reactivity:
human

GenBank:
BC012061

Immunogen Catalog Number:
Ag9202

Conjugate:
Unconjugated

Full name:
poly(rC) binding protein 3

Gene ID:
54039

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01894-1 targets PCBP3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PCBP3 Recombinant antibody, PBS Only (Capture) 85180-2-PBS (242791E7). 100 µg. Concentration 1 mg/mL.

Detection antibody: PCBP3 Recombinant antibody, PBS Only (Capture/Detector) 85180-3-PBS (242791F11). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01894-1, PCBP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85180-2-PBS. Detection antibody: 85180-3-PBS. Standard: Ag9202. Range: 0.781-100 ng/mL.