

For Research Use Only



KAP1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP02033-1

Capture Antibody Information

Catalog Number:
85693-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
243043C1
Reactivity:
human
GenBank:
BC004978
Immunogen Catalog Number:
Ag7350

Conjugate:
Unconjugated
Full name:
tripartite motif-containing 28
Gene ID:
10155

Detection Antibody Information

Catalog Number:
85693-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
243043A7
Reactivity:
human
GenBank:
BC004978
Immunogen Catalog Number:
Ag7350

Conjugate:
Unconjugated
Full name:
tripartite motif-containing 28
Gene ID:
10155

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02033-1 targets KAP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KAP1 Recombinant antibody, PBS Only (Capture) 85693-3-PBS (243043C1). 100 µg. Concentration 1 mg/mL.

Detection antibody: KAP1 Recombinant antibody, PBS Only (Detector) 85693-2-PBS (243043A7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02033-1, KAP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85693-3-PBS. Detection antibody: 85693-2-PBS. Standard: Ag7350. Range: 0.625-80 ng/mL