

For Research Use Only

# IL6 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50032-1

## Capture Antibody Information

**Catalog Number:**  
68659-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
6C12C4  
**Reactivity:**  
Human  
**GenBank:**  
Immunogen Catalog Number:  
HZ-1019

**Conjugate:**  
Unconjugated  
**Full name:**  
interleukin 6 (interferon, beta 2)  
**Gene ID:**  
3569

## Detection Antibody Information

**Catalog Number:**  
68659-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
6C4B7  
**Reactivity:**  
Human  
**GenBank:**  
Immunogen Catalog Number:  
HZ-1019

**Conjugate:**  
Unconjugated  
**Full name:**  
interleukin 6 (interferon, beta 2)  
**Gene ID:**  
3569

## Applications

**Tested Applications:**  
Sandwich ELISA

**Range:**  
7.8-500 pg/mL (Sandwich ELISA)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50032-1 targets IL6 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: IL-6 Monoclonal antibody, PBS Only (Capture) 68659-1-PBS (6C12C4). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL-6 Monoclonal antibody, PBS Only (Detector) 68659-2-PBS (6C4B7). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP50032-1, IL6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68659-1-PBS. Detection antibody: 68659-2-PBS. Standard: HZ-1019. Range: 7.8-500 pg/mL